

Implant compatibility chart

Titanium and gold-shaded titanium

Internal connection interfaces		
Manufacturer	Implant type and diameter	Abutment screw type
Astra Tech	OsseoSpeed™ 3.0 (Yellow) OsseoSpeed™ 3.5/4.0 (Aqua) OsseoSpeed™ 4.5/5.0 (Lilac)	0.050" Hex Drive
BioHorizons	Internal 3.5, 4.5, 5.7	1.25 mm Hex Drive
Biomet 3i	Certain™ MicroMiniplant™ 3.25/3.4 Certain™ 3.75/4.0, 5.0, 6.0, XP 4/5, XP 5/6	0.048" Hex Drive
Nobel Biocare	NobelReplace™ NP (Narrow Platform 3.5) NobelReplace™ RP (Regular Platform 4.3) NobelReplace™ WP (Wide Platform 5.0, 6.0)	Nobel Biocare Unigrip™
Straumann	Standard/Standard Plus Implant RN (Regular Neck) 4.8	Straumann Drive Compatible
Zimmer Dental	Tapered Screw-Vent™ 3.5, 4.5, 5.7 Screw-Vent™ 3.5, 4.5 SwissPlus™ 4.8	0.050" Hex Drive or 1.25 mm Hex Drive

External hex interfaces		
Manufacturer	Implant type and diameter	Abutment screw type
BioHorizons	External 3.5, 4.0, 5.0	0.050" Hex Drive
Biomet 3i	MicroMiniplant™ 3.25/3.4 Miniplant™ 4.1/3.25 Standard 3.75, 4.0 Wide 5.0, 6.0 XP 3/4, 4/5, 5/6	0.048" Hex Drive
Nobel Biocare	Brånemark System™ NP (Narrow Platform 3.3) Brånemark System™ RP (Regular Platform 3.75, 4.0) Brånemark System™ WP (Wide Platform 5.0, 5.5)	Nobel Biocare Unigrip™

All orders for Atlantis™ abutments are shipped with the corresponding abutment screw needed at no extra charge.

Note: The implant manufacturer's recommended torque should always be used when placing the abutment screw.

Turn around time

Units per case	Turn around time (in-house business days)	
	Titanium	GoldHue™
1-3	4*	10
4-5	5	11
6+	6	12

Note: Turn around time is from door-to-door, the day it is shipped from the customer to the day it returns to the customer.

To meet the turn around times the shipping from the customer has to be next-day noon delivery, shipping back to the customer will always be next day delivery.

Turn around times do not include holidays and weekends.

GoldHue™ has a titanium nitride coating. It has a longer turn around time due to the nature of the coating process.

*E.g. if 1-3 units titanium is shipped from the customer on Monday, the customer will receive the order on Friday the same week.

Internal connection interfaces		
Manufacturer	Implant type and diameter	Abutment screw type
Astra Tech	OsseoSpeed™ 3.5/4.0 (Aqua) OsseoSpeed™ 4.5/5.0 (Lilac)	0.050" Hex Drive
BioHorizons	Internal 3.5, 4.5, 5.7	1.25 mm Hex Drive
Biomet 3i	Certain™ MicroMiniplant™ 3.25/3.4 Certain™ 3.75/4.0, 5.0, 6.0, XP 4/5, XP 5/6	0.048" Hex Drive
Nobel Biocare	NobelReplace™ RP (Regular Platform 4.3) NobelReplace™ WP (Wide Platform 5.0, 6.0)	Nobel Biocare Unigrip™
Zimmer Dental	Tapered Screw-Vent™ 3.5, 4.5, 5.7 Screw-Vent™ 3.5, 4.5	0.050" Hex Drive or 1.25 mm Hex Drive

External hex interfaces		
Manufacturer	Implant type and diameter	Abutment screw type
BioHorizons	External 4.0, 5.0	0.050" Hex Drive
Biomet 3i	Miniplant™ 4.1/3.25 Standard 3.75, 4.0 Wide 5.0, 6.0 XP 3/4, 4/5, 5/6	0.048" Hex Drive
Nobel Biocare	Brånemark System™ RP (Regular Platform 3.75, 4.0) Brånemark System™ WP (Wide Platform 5.0, 5.5)	Nobel Biocare Unigrip™

All orders for Atlantis™ abutments are shipped with the corresponding abutment screw needed at no extra charge.

Note: The implant manufacturer's recommended torque should always be used when placing the abutment screw.

Turn around time

Units per case	Turn around time (in-house business days)
	Zirconia
1-3	7
4-5	8
6+	9

Note: Turn around time is from door-to-door, the day it is shipped from the customer to the day it returns to the customer.

To meet the turn around times the shipping from the customer has to be next-day noon delivery, shipping back to the customer will always be next day delivery.

Turn around times do not include holidays and weekends.

All trademarks are the property of their respective owners.

